REF. 2138 289 89 K

| | | |
|-----------|--|-----------|
| FR | NOTICE D'UTILISATION - SONDES BICANALICULONASALES | 2 |
| EN | INSTRUCTION FOR USE - BICANALICULAR INTUBATION SET | 3 |
| DE | GEBRAUCHSANWEISUNG - BIKANALIKULÄRES INTUBATIONSET | 4 |
| IT | ISTRUZIONI PER L'USO - SONDE NASALI BICANALICOLARI | 5 |
| ES | INSTRUCCIONES DE USO - SONDAS BICANALICULONASALES | 6 |
| PT | INSTRUÇÕES DE USO - SONDAS BICANALICULONASAIS | 7 |
| NL | GEBRUIKSAANWIJZING - BICANALICULAIRE NEUSSONDES | 8 |
| DA | BRUGSVEJLEDNING - BIKANALIKULÆRT INTUBATIONSSÆT | 9 |
| FI | KÄYTTÖOHJEET - KAKSIKANAVAINEN SONDEERAUSSETTI | 10 |
| SE | BRUKSANVISNING - BIKANALIKULÄRT INTUBATIONSSET | 11 |
| TR | KULLANIM - BIKANALİKÜLER ENTÜBASYON SETİ | 12 |
| ZH | 使用说明 - 泪道引流管 (BIKA 双管型) | 13 |
| EL | ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΣΕΤ ΔΙΣΩΛΗΝΑΡΙΑΚΗΣ ΔΙΑΣΩΛΗΝΩΣΗΣ | 15 |
| PL | INSTRUKCJA UŻYTKOWANIA – ZESTAW DO INTUBACJI DWUKANALIKOWEJ | 16 |
| SL | NAVODILA ZA UPORABO – KOMPLET ZA BIKANALIKULARNO INTUBACIJO | 17 |
| SK | NÁVOD NA POUŽITIE - SÚPRAVA NA BIKANALIKULÁRNU INTUBÁCIU | 18 |

Symbols..... 19

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EN

Instructions for use
Date of revision of the instructions for use : 04/2016

BICANALICULAR INTUBATION SET

Description :

BICANALICULAR INTUBATION SETS are composed of a silicone tube connected by each of its extremities to a malleable steel guide. The silicone tube may be coated with polyvinylpyrrolidone (PVP) to improve its wettability. There are different models of BICANALICULAR INTUBATION SETS depending on:

- The tube to metal connection
- The size of the silicone tube
- The size of the metal guides

Packaging :

BICANALICULAR INTUBATION SETS are supplied sterile. They are sterilized by ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green. BICANALICULAR INTUBATION SETS are sterilized in their final double-blister or double-pouch packaging to make them easier to handle in aseptic conditions.

Indications :

Nasal bicanalicular intubation is indicated in treatments of epiphora in infants or adults, particularly in cases of:

- Canalicular pathologies (stenosis, obstruction, lacerations),
- Dacryocystorhinostomy (conventional or laser),
- Imperforation of the nasolacrimal duct in the infants.

Mode of action :

Nasal bicanalicular intubations are used for intubation of the lacrimal ducts. The silicone tube acts as a conformer and enables drainage of tears by capillarity. In cases of canalicular lacerations, the silicone tube guides wound healing and prevents the onset of synechia.

Warning and contraindications :

The use of BICANALICULAR INTUBATION SETS is contra-indicated in cases of epiphora that do not have a nasolacrimal origin. The silicone / steel guide junction is a fragile area. It is recommended that the steel guide be held in front of the junction while handling the probe.

Adverse effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of the BICANALICULAR INTUBATION SETS include, but are not limited to the following :

Intraoperative complications :

- false passages
- separation of the silicone from the steel guide

Post-operative complications :

- conjunctival or nasal pruritus
- nasal or caruncular irritation
- reversible shrinkage of the palpebral fissure
- exteriorization or loss of the probe
- stricturotomy
- synechia of the nasal mucosa
- canaliculitis
- induced mucocele

The application of eye drops containing methylene blue can cause staining of the tube. This has no effect on the patient's tolerability of the probe. Unexpected side effects and complications related to the BICANALICULAR INTUBATION SETS must be reported to FCI.

Recommendations for use :

It is recommended that the patient be shown the interpalpebral loop after the intervention.

Precautions for use :

BICANALICULAR INTUBATION SETS must be removed from their packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. BICANALICULAR INTUBATION SETS are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. They should be stored at room temperature and must not be used after the expiration date shown on the package.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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