REF. 2138 289 93 K

FR	NOTICE D'UTILISATION - CLOUS-TROUS® PVP	2
EN	INSTRUCTION FOR USE - PVP PERFORATED PLUGS	3
DE	GEBRAUCHSANWEISUNG - PERFORIERTER PVP PLUG	4
IT	ISTRUZIONI PER L'USO - TAPPI-FORI® PVP	5
ES	INSTRUCCIONES DE USO - TAPONES-PERFORADOS® PVP	6
PT	INSTRUÇÕES DE USO - ANILHAS PERFURADAS® PVP	7
NL	GEBRUIKSAANWIJZING - PERFORATUS PLUG® PVP	8
DA	BRUGSVEJLEDNING - PVP PERFOREREDE PROPPER	9
FI	KÄYTTÖOHJEET - REIÄLLISET PVP-TULPAT	10
SE	BRUKSANVISNING - PVP-PERFORERAD PLUGG	11
TR	KULLANIM - PVP DELİKLİ TİKAÇLAR	12
ZH	使用说明 - 泪道引流管 (Perforated PlugS 穿孔塞型)	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΔΙΑΤΡΗΤΑ ΠΩΜΑΤΑ ΜΕ ΕΠΕΝΔΥΣΗ PVP	15
PL	INSTRUKCJA UŻYTKOWANIA – ZATY CZKI PERFOROWANE PVP	16
SL	NAVODILA ZA UPORABO – PERFORIRANI ČEPI PVP	17
SK	NÁVOD NA POUŽITIE - PERFOROVANÉ ZÁTKY PVP	18

Symbols 19

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
PVP PERFORATED PLUGS

Description and packaging:

PVP PERFORATED PLUGS are small molded pieces of silicone that consist of:

- a slanted collarette
- a body
- a bulb.

PVP PERFORATED PLUGS are coated with polyvinylpyrrolidone (PVP) to improve their wettability and thus limit risk of occlusion.

PVP PERFORATED PLUGS are supplied sterile. They are sterilized in their final packaging by gamma irradiation (minimal dose: 25 kGy). A red indicator shows that the product followed a sterilization cycle validated by FCI.

The PVP PERFORATED PLUGS are individually packaged and preloaded on a disposable PLUG INSERTER ("PRE-LOADED" PERFORATED PLUGS), equipped with dilator and packed in single-blister packaging.

Indications:

PVP PERFORATED PLUGS are used for draining tears in the treatment of epiphora caused by a punctal stenosis, while the rest of the lacrimal duct is normal.

Mode of action:

The PERFORATED PLUG acts as a conformer which gradually dilates the punctum. This dilation is durable after implant's removal.

The three parts of the PVP PERFORATED PLUG perform the following functions:

The collarette covers the punctal opening and keeps the plug from migrating. It is thin and slanted to fit the eyelid margin thus limiting the risk of corneal irritation, The bulb is fixed in the vertical part of the canaliculus and prevents spontaneous expulsion of the PLUG,

The body links the collarette to the bulb and includes an open canal to enable drainage of tears.

Contraindications:

Implantation of PVP PERFORATED PLUG is contra-indicated in the following cases:

- Impermeability of the lacrimal ducts,
- Abnormalities of the eyelids with absence of nocturnal palpebral occlusion,
- Senile conjunctival hyperlaxity covering the inferior punctum,
- Canalicular stenosis,
- Ruptured punctum.

Adverse side effects:

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of PVP PERFORATED PLUGS include, but are not limited to the following :

Complications occurring during insertion of PVP PERFORATED PLUGS :

- Pain due to dilation of the punctum,
- Rupture of the meatus.

Postoperative complications :


- Temporary pruritus lasting for approximately a week,
- Temporary dryness,
- Chronic irritation,
- Expulsion of the PLUG,
- Pyogenic granuloma,
- Canalicular stenosis.

Unexpected side effects and complications related to PVP PERFORATED PLUG must be reported to FCI.

Instructions for use:

Insertion:

- topical anesthesia of the eyelid,
- careful dilation of the meatus with the dilator of the PLUG INSERTER,
- insertion of the PLUG using a PLUG INSERTER.
- Patient should be relaxed, seated or supine
- Proper lighting is essential for good visibility
- Although the entire procedure takes only a few minutes, it should not be rushed
- Generously dilate the punctum
- Grasp eyelid with one hand to stabilize it horizontally and expose the punctal opening
- Before inserting, orient the plug so that the collarette slants up nasally (in inferior punctum) and down nasally (in superior punctum).
- Introduce plug vertically, in a gentle rocking and progressive manner. The inserter will retreat slightly at the moment the bulb portion of the plug goes through the punctal opening.
- Once properly positioned, release plug by squeezing the instrument handles together
- Examine patient with slit lamp
- If necessary, rotate plug so its collarette lies flush against eyelid margin
- To facilitate insertion, it may be helpful to lubricate the plug with saline, artificial tears, or a similar solution.
- The PLUG should usually be left in place for 2 to 3 months (this period is long enough for the treatment to be effective and short enough to avoid complications).


Removal : Removal is painless. Use blunt forceps to grasp PLUG below the collarette and gently pull it out from punctal opening.

Warnings and precautions for use:

PUNCTAL PLUGS must be removed from their blister and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. PUNCTAL PLUGS are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not red as it might endanger the patient. They should be stored at room temperature and must not be used after the expiration date shown on the package.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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