**REF. 2138 289 94 L**

<b>FR</b>	<b>NOTICE D'UTILISATION - BOUCHONS MEATIQUES</b>	<b>2</b>
<b>EN</b>	<b>INSTRUCTION FOR USE - PUNCTAL PLUGS</b>	<b>3</b>
<b>DE</b>	<b>GEBRAUCHSANWEISUNG - PUNCTUM-PLUG</b>	<b>4</b>
<b>IT</b>	<b>ISTRUZIONI PER L'USO - TAPPI MEATICI</b>	<b>5</b>
<b>ES</b>	<b>INSTRUCCIONES DE USO - TAPONES MEÁTICOS</b>	<b>6</b>
<b>PT</b>	<b>INSTRUÇÕES DE USO - TAMPÕES MEÁTICOS</b>	<b>7</b>
<b>NL</b>	<b>GEbruIKSAANWIJZING - DOORGANGSPLUG</b>	<b>8</b>
<b>DA</b>	<b>BRUGSVEJLEDNING - PUNCTALPROPPER</b>	<b>9</b>
<b>FI</b>	<b>KÄYTTÖOHJEET - PUNCTUM-TULPAT</b>	<b>10</b>
<b>SE</b>	<b>BRUKSANVISNING - PUNCTUMPLUGG</b>	<b>11</b>
<b>TR</b>	<b>KULLANIM - PUNKTUM TIKAÇLARI</b>	<b>12</b>
<b>ZH</b>	<b>使用说明 - 泪点塞 (预装式)</b>	<b>13</b>
<b>EL</b>	<b>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΠΩΜΑΤΑ ΔΑΚΡΥΪΚΟΥ ΣΗΜΕΙΟΥ</b>	<b>15</b>
<b>PL</b>	<b>INSTRUKCJA UŻYTKOWANIA – ZATYCZKI PUNKTÓW ŁZOWYCH</b>	<b>16</b>
<b>SL</b>	<b>NAVODILA ZA UPORABO – ČEPI ZA SOLZNO ODPRTINO VEKE</b>	<b>17</b>
<b>SK</b>	<b>NÁVOD NA POUŽITIE - ZÁTKY SLZNÝCH BODOV</b>	<b>18</b>

**Symbols..... 19**

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
## PUNCTAL PLUGS

### Description and packaging:

PUNCTAL PLUGS are small molded pieces of silicone that consist of:

- a slanted collarette
- a body
- a bulb.

PUNCTAL PLUGS are available in different sizes.

PUNCTAL PLUGS are supplied sterile. They are sterilized in their final packaging by gamma irradiation (minimal dose: 25 kGy). A red indicator shows that the product followed a sterilization cycle validated by FCI.

The PUNCTAL PLUGS are individually packed mounted on a disposable PLUG INSERTER (« PRE-LOADED » PUNCTAL PLUGS) equipped with dilator and packed in a single packaging.

### Indications :

PUNCTAL PLUGS are designed to close the lacrimal puncta in the treatment of dry eye. The closure is more effective when both puncta of a single side are closed.

### Mode of action:

The three parts of the PUNCTAL PLUG fulfill the following functions:

The collarette covers the punctal opening and keeps the plug from migrating. It is thin and slanted on the eyelid margin thus limiting the risk of corneal laceration. The bulb is fixed in the vertical part of the canaliculus and prevents spontaneous expulsion of the PUNCTAL PLUG. The body links the collarette to the bulb and includes a closed canal for plug placement.

### Contraindications:

Implantation of PUNCTAL PLUGS is contraindicated in the following cases:

- Impermeability of the lacrimal ducts,
- Abnormalities of the eyelids with absence of nocturnal palpebral occlusion,
- Senile conjunctival hyperlaxity covering the inferior meatus,
- A positive Bengal Rose test beyond the area of the palpebral fissure,
- Rupture of the meatus.

### Adverse side effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of the PUNCTAL PLUGS include, but are not limited to the following:

#### Complications occurring during insertion of PUNCTAL PLUGS:

- Pain due to dilation of the punctum,
- Rupture of the meatus.

#### Complications occurring when the PUNCTAL PLUG is in place:

- Temporary pruritus lasting for approximately a week,
- Epiphora,
- Chronic irritation,
- Expulsion of the punctal plug,
- Canalicular stenosis occurring after loss of the punctal plug,
- Pyogenic granuloma.

Unexpected side effects and complications related to PUNCTAL PLUGS must be reported to FCI.


### Instructions for use :

It is recommended to use the PUNCTAL PLUG GAUGES (S2-3061, S2-3062, S2-3063 and S2-3064) to determine the adequate size of the plug.

#### Insertion:

A topical anesthesia of the eyelid is recommended, excepted for the EXTRA-SMALL PLUGS, carefully dilate the meatus with the dilator of the PLUG INSERTER, insert the PUNCTAL PLUG using a PLUG INSERTER

- Patient should be relaxed, seated or supine
- Proper lighting is essential for good visibility
- Although the entire procedure takes only a few minutes, it should not be rushed
- Generously dilated punctum
- Grasp eyelid with one hand to stabilize it horizontally and expose the punctal opening
- Before inserting, orient the plug so that the collarette slants up nasally (in inferior punctum) and down nasally (in superior punctum).
- Likeness of the plug stamped on the insertion instrument shows the direction of the slant for proper orientation
- Introduce plug vertically, in a gentle rocking and progressive manner. The inserter will retreat slightly at the moment the bulb portion of the plug goes through the punctal opening.
- Once properly positioned, release plug by squeezing the instrument handles together
- Examine patient with slit lamp
- If necessary, rotate plug so its collarette lies flush against eyelid margin
- To facilitate insertion, it may be helpful to lubricate the plug with saline, artificial tears, or a similar solution.

The plug should be left in place for 2 or 3 months (this period is long enough for the treatment to be effective and short enough to avoid complications).


#### Removal:

Removal is painless. Use blunt forceps to grasp plug below the collarette and gently pull it out from punctal opening.

### Warnings and precautions for use:

PUNCTAL PLUGS must be removed from their packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. PUNCTAL PLUGS are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not red as it might endanger the patient. They should be stored at room temperature and must not be used after the expiration date shown on the package.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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