REF. 2138 289 99 I

FR	NOTICE D'UTILISATION - BILLES EN SILICONE	2
EN	INSTRUCTION FOR USE - SILICONE ORBITAL IMPLANTS	3
DE	GEBRAUCHSANWEISUNG - ORBITAIMPLANTAT AUS SILIKON	4
IT	ISTRUZIONI PER L'USO - SFERE DI SILICONE	5
ES	INSTRUCCIONES DE USO - BOLAS DE SILICONA	6
PT	INSTRUÇÕES DE USO - ESFERAS DE SILICONE	7
NL	GEbruIKSAANWIJZING - ORBITALE SILICONENIMPLANTAAT	8
DA	BRUGSVEJLEDNING - ORBITALIMPLANTATER I SILIKONE	9
FI	KÄYTTÖOHJEET - SILIKONISET SILMÄKUOPPAIMPLANTIT	10
SE	BRUKSANVISNING - ORBITAIMPLANTAT AV SILIKON	11
TR	KULLANIM - SİLİKON ORBİTAL İMPLANTLAR	12
ZH	使用说明 - 硅胶型眼眶植入物	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΕΝΔΟΚΟΓΧΙΚΑ ΕΜΦΥΤΕΥΜΑΤΑ ΣΙΛΙΚΟΝΗΣ	14
PL	INSTRUKCJA UŻYTKOWANIA – IMPLANTY OCZODOŁOWE SILIKONOWE	15
SL	NAVODILA ZA UPORABO – SILIKONSKI ORBITALNI VSADKI	16
SK	NÁVOD NA POUŽITIE - SILIKÓNOVÝ ORBITÁLNY IMPLANTÁT	17

Symbols..... 18

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
SILICONE ORBITAL IMPLANTS

Description and packaging:

ENUCLEATION / EVISCERATION IMPLANTS are orbital implants designed to fill the orbital cavity following enucleation or evisceration.

ENUCLEATION / EVISCERATION IMPLANTS are made of silicone and are available in the following diameters: 12, 13, 14, 15, 16, 17, 18, 19, 20 and 22 mm (S6.1012, S6.1013, S6.1014, S6.1015, S6.1016, S6.1017, S6.1018, S6.1019, S6.1020, S6.1022).

ENUCLEATION / EVISCERATION IMPLANTS are supplied sterile and sterilized in their final double-blister packaging, making them easier to handle in aseptic conditions. ENUCLEATION / EVISCERATION IMPLANTS are sterilized by gamma irradiation (minimal dose: 25 kGy). A red indicator shows that the product followed a sterilization cycle validated by FCI.

Indications:

ENUCLEATION / EVISCERATION IMPLANTS are orbital implants designed to fill the orbital cavity following enucleation, evisceration or secondary implantation.

Complications:

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of the ENUCLEATION / EVISCERATION IMPLANTS include, but are not limited to the following :

- Conjunctival dehiscence
- Localized infection or irritation
- Exposure
- Discomfort
- Ptosis
- Laxity of eyelids
- Enophthalmos
- Hypophthalmos
- Implant expulsion

Unexpected side effects and complications related to the ENUCLEATION / EVISCERATION IMPLANTS must be reported to FCI.

Instructions for use :

The choice of the diameter is an important parameter for the final result. To avoid dehiscence around the incision, the diameter of ENUCLEATION/EVISCERATION IMPLANTS has to be adapted to the patient's anatomy.

Sizers S6.3060 are used to determine the appropriate implant diameter. The proper diameter should allow the wound to be closed without tension.

Warnings and precautions for use:

Remove the ENUCLEATION/EVISCERATION IMPLANTS from their packaging in aseptic conditions. Before use, verify that the individual pack preserving the product's sterility should be checked to make sure it is intact.

ENUCLEATION/EVISCERATION IMPLANTS are single-use products that should not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not red, as it could endanger the patient. They should be stored at room temperature and should not be used after the expiration date indicated on the package.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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