REF. 2138 290 00 G

FR	NOTICE D'UTILISATION - BILLES HA®	2
EN	INSTRUCTION FOR USE - HA ORBITAL IMPLANTS	3
DE	GEBRAUCHSANWEISUNG - HA-ORBITAIMPLANTAT	4
IT	ISTRUZIONI PER L'USO - SFERE HA®	5
ES	INSTRUCCIONES DE USO - BOLAS HA®	6
PT	INSTRUÇÕES DE USO - ESFERAS HA®	7
NL	GEBRUIKSAANWIJZING - HA® IMPLANTATEN	8
DA	BRUGSVEJLEDNING - HA ORBITALIMPLANTATER	9
FI	KÄYTTÖOHJEET - HA SILMÄKUOPPAIMPLANTIT	10
SE	BRUKSANVISNING - HA-ORBITAIMPLANTAT	11
TR	KULLANIM - HA ORBİTAL İMPLANTLAR	12
ZH	使用说明 - HA 眼眶植入物	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΕΝΔΟΚΟΓΧΙΚΑ ΕΜΦΥΤΕΥΜΑΤΑ HA	14
PL	INSTRUKCJA UŻYTKOWANIA – IMPLANTY OCZODOŁOWE HA	15
SL	NAVODILA ZA UPORABO – HA ORBITALNI VSADKI	16
SK	NÁVOD NA POUŽITIE - HA ORBITÁLNE IMPLANTÁTY	17

Symbols 18

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
HA ORBITAL IMPLANTS

Description and packaging:

HA ORBITAL IMPLANTS are intra-orbital synthetic hydroxyapatite ($\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$) implants with 95 % guaranteed purity. Their porosity is range from 300 to 500 μm . They are available in several diameters from 10 mm to 22 mm (S6.2010, S6.2012, S6.2014, S6.2016, S6.2018, S6.2020, S6.2022).

HA ORBITAL IMPLANT

The HA ORBITAL IMPLANTS are supplied sterile. They are sterilized using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. HA ORBITAL IMPLANTS are sterilized in their final double packaging, to make them easier to handle in aseptic conditions.

Indications:

The HA ORBITAL IMPLANTS are designed to fill in the orbital cavity following enucleation, evisceration or during secondary implantation.

Mode of action:

The porosity of the hydroxyapatite facilitates its colonization by fibrovascular tissue, which offers the following advantages:

- reduced risk of infection
- reduced risk of implant extrusion

Contraindications:

Care should be given to irradiated orbital cavities where the capacity for vascularization and colonization can be reduced (especially in cases of irradiation of the orbital cavity and in diabetic patients).

Complications:

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of the HA ORBITAL IMPLANTS include, but are not limited to the following :

- | | |
|-------------------------------------|---------------------|
| - Conjunctival dehiscence | - Laxity of eyelids |
| - Localized infection or irritation | - Enophthalmos |
| - Exposure | - Hypophthalmos |
| - Discomfort | |
| - Ptosis | |

Unexpected side effects and complications related to the HA ORBITAL IMPLANTS must be reported to FCI.

Intructions for use :

The choice of the diameter is an important parameter for the final result. Given the retractile effect of orbital tissues, and to avoid dehiscence around the scar, the HA ORBITAL IMPLANT diameter has to be adapted to the patient's anatomy.

The use of sizers S6.3060 is advisable. The proper diameter should allow the wound to be closed without tension. The choice of the diameter must take into account for the thickness of any material used to wrap the implant.

The use of sphere introducer S6.3050 is recommended for the install of the HA ORBITAL IMPLANT.

The enucleation or the evisceration technique is conventional

HA ORBITAL IMPLANTS can be implanted uncovered, particularly in cases of evisceration.

In cases of evisceration, it is usually not advisable to retain the cornea, particularly if there is history of corneal injury.

In cases of enucleation :

Or cover the HA ORBITAL IMPLANTS with either absorbable mesh or autogenous material (sclera, fascialata, temporal aponeurosis)

Warnings and precautions for use:

Remove the HA ORBITAL IMPLANTS from their packaging in aseptic conditions. Before use, verify that the individual pack preserving the product's sterility should be checked to make sure it is intact.

HA ORBITAL IMPLANTS are single-use products that should not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green, as it could endanger the patient. They should be stored at room temperature and should not be used after the expiration date indicated on the package.

The drilling of the sphere is not recommended as the process may damage it.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 μg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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