**REF. 2138 291 84 F**

FR	NOTICE D'UTILISATION – PTOSE-UP®	2
EN	INSTRUCTION FOR USE – PTOSE-UP®	3
DE	GEBRAUCHSANWEISUNG – PTOSE-UP®	4
IT	ISTRUZIONI PER L'USO – PTOSE-UP®	5
ES	INSTRUCCIONES DE USO – PTOSE-UP®	6
PT	INSTRUÇÕES DE USO – PTOSE-UP®	7
NL	GEBRUIKSAANWIJZING – PTOSE-UP®	8
DA	BRUGSVEJLEDNING – PTOSE-UP®	9
FI	KÄYTTÖOHJEET – PTOSE-UP®	10
SE	BRUKSANVISNING – PTOSE-UP®	11
TR	KULLANIM – PTOSE-UP®	12
ZH	使用说明 – 下垂上拉带	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ – PTOSE-UP®	15
PL	INSTRUKCJA UŻYTKOWANIA – PTOSE-UP®	16
SL	NAVODILA ZA UPORABO – PTOSE-UP®	17
SK	NÁVOD NA POUŽITIE – PTOSE-UP®	18

**Figures.....19**

**Symbols.....20**

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## PTOSE-UP®

### Description :

The PTOSE-UP® is made of ePTFE (expanded polytetrafluoroethylene) and is a ready-to-use pre-cut band available in two widths: 3 mm or 2 mm. Each extremity of the band is fitted with a needle for its installation by tunneling technique.

### Packaging :

The PTOSE-UP® is sterilized in its final packaging with ethylene oxide: a green indicator shows that the product followed a sterilization cycle validated by FCI S.A.S. Do not use the product if the indicator is not green. The PTOSE-UP® is sterilized in its final double packaging to make it easier to handle in aseptic conditions.

### Indications :

The PTOSE-UP® is used for the surgical treatment of ptosis requiring brow suspension, where a Levator palpebrae is unusable following a functional deficiency or anatomic destruction, particularly in the following cases:

- certain major congenital types of ptosis
- certain types of ptosis associated with blepharophimosis syndrome
- certain acquired myogenic and neurogenic forms of ptosis
- certain cases of post-traumatic ptosis
- certain forms of ptosis with innervational disorders (e. g. Marcus Gunn syndrome)
- certain forms of III paralysis

### Mode of action:

In cases where the Levator palpebrae is unusable, the surgical treatment of ptosis requires brow suspension, with the band being placed between the tarsal plate and the occipito-frontalis muscle. The PTOSE-UP® is sutured between the upper edge of the tarsus and the levator muscle with resorbable 6.0 suture. Using the ePTFE PTOSE-UP® rather than autogenous or homologous (fascia lata, temporal aponeurosis,...) offers the following particular advantages:

- reduced operating and hospitalization time
- a single operative site
- elimination of the inherent morbidity linked to harvesting autogenous material
- biocompatibility and biointegration of the ePTFE PTOSE-UP®

### Figure 2:

- 1 – Dissection of a flap of pre-tarsal orbicularis (below the clamps) which is used to protect the biomaterial from the cutaneous incision.
- 2 – Closure of the eyelid incision.
- 3 – Suturing the ePTFE to the frontalis muscle.
- 4 – Suturing the supra-eyebrow incision at two levels (the upper offset from the line of the deeper suture, which is 4 – 5 mm lower).

### Warning and contra-indication:

In cases of bilateral ptosis, where there are risks of corneal trophic disorders induced by secondary lagophthalmos following suspension (pre-existing corneal hypoaesthesia, absence or limitation of the Charles Bell sign, certain major myogenic or neurogenic ptoses), it is recommended that the patient be operated upon in 2 stages with a 30-45 days interval in between.

The PTOSE-UP® shall not be manipulated with toothed forceps as the instrument could damage the material. Do not use powdered surgical gloves as the direct contact of polymer with powder might obstruct the permeability of the porous material and modify its properties.

The surgical procedure requires very strict aseptic conditions (routine use of polyvidone iodide solution) and the greatest possible protection of the biomaterial from the cutaneous incision (where there is a risk of its extrusion before it is colonized by cells) by a layer of high quality, well vascularized tissue.

### Adverse side effects:

As in any type of surgery, there are risks linked to the material and to developments in the initial pathology. Potential complications associated with surgery for ptosis using the PTOSE-UP® include, but are not limited to the following :

- corneal ulcers by post-operative lagophthalmos, particularly in patients with myogenic or neurogenic ptosis
- delayed scar formation, granulomas, or extrusion of the biomaterial at the area and level of the incisions
- under- or over-corrections

Side effects or complications reasonably thought to be caused by the PTOSE-UP® must be reported to FCI S.A.S.

### Recommendations for use:

The PTOSE-UP® should be removed from its package just before the procedure.

A Frost suture is recommended to protect the cornea at the end of the surgery.

Disinfecting the operative site and soaking the PTOSE-UP® during the intervention with a solution of iodized polyvidone are recommended.

### Warnings :

The PTOSE-UP® must be removed from its packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. The PTOSE-UP® is a single-use product and must not be resterilized. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety. It should be stored at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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