REF. 2138 295 57 D

FR	NOTICE D'UTILISATION - POSE-CLOU DILATATEUR 0.3 A USAGE UNIQUE	2
EN	INSTRUCTION FOR USE - SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER	3
DE	GEBRAUCHSANWEISUNG - PLUG-INSERTER UND PUNCTUM-DILATOR 0.3 ZUM EINMALGEBRAUCH	4
IT	ISTRUZIONI PER L'USO - POSA-TAPPO DILATATORE 0.3 USA E GETTA	5
ES	INSTRUCCIONES DE USO - APLICADOR DILATADOR 0.3 PARA USO ÚNICO	6
PT	INSTRUÇÕES DE USO - APLICADOR-PREGO DILATADOR 0.3 PARA USO ÚNICO	7
NL	GEBRUIKSAANWIJZING - APPLICATOR EN DILATATOR 0,3 VOOR EENMALIG GEBRUIK	8
DA	BRUGSVEJLEDNING - PUNKTALDILATOR 0,3 TIL ENGANGSBRUG MED IND-FØRINGSANORDNING	9
FI	KÄYTTÖOHJEET - KERTAKÄYTTÖINEN 0.3 PUNKTUMLAAJENTAJA JA TULPPA-ASETIN	10
SE	BRUKSANVISNING - PUNCTUM-DILATOR OCH PLUG-INSERTER 0.3 FÖR ENGÅNGSBRUK	11
TR	KULLANIM - TEK KULLANIMLIK 0.3 PUNCTUM DİLATORÜ VE TIKAÇ YERLEŞTİRİCİSİ	12
ZH	使用说明 - 一次性 0.3 泪点扩张器及插塞器	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΔΙΑΣΤΟΛΕΑΣ ΔΑΚΡΥΪΚΟΥ ΣΗΜΕΙΟΥ ΚΑΙ ΕΙΣΑΓΩΓΕΑΣ ΠΩΜΑΤΟΣ 0,3 ΜΙΑΣ ΧΡΗΣΗΣ	14
PL	INSTRUKCJA UŻYTKOWANIA – JEDNORAZOWY ROZSZERZACZ PUNKTU ŁZOWEGO 0.3 I APLIKATOR ZATYCZKI	15
SL	NAVODILA ZA UPORABO – DILATATOR SOLZNE LUKNJICE 0,3 IN VSTAVLJALEC ČEPOV ZA ENKRATNO UPORABO	16
SK	NÁVOD NA POUŽITIE - JEDNORAZOVÝ 0,3 DILATÁTOR SLZNÉHO BODU A ZAVÁDZAČ ZÁTKY	17

Symbols..... 18

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
SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER

Description and packaging:

The SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER is made of polypropylene with

- a metallic tip 0.3 mm in diameter, the same diameter as the channel inside the fixation head (plug) of the MONOKA® monocanalicular intubation sets
- a mechanism that allows a progressive release by squeezing with the thumb and index finger
- a dilator


The SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER is packaged in a double pouch and supplied sterile. It is sterilized in its final packaging by gamma irradiation (minimal dose: 25 kGy). A red indicator shows that the product followed a sterilization cycle validated by FCI.

Indications :

The SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER facilitates the insertion of the fixation head (plug) of MONOKA® monocanalicular intubation sets.

Instructions for use :

- After retrieval of the silicone tube of the monocanalicular intubation set in the lacrimal duct, pull on it until the plug at the opposite end reaches the punctum
- Hold the SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER between thumb and index finger at the level of the mechanism without squeezing.
- Insert the metallic tip into the hollow body of the plug that links the collarette to the bulb.
- Push the plug into the punctal opening
- Squeeze the instrument handles together to retract the insertion tip leaving the plug in position in the punctum.

Warnings and precautions for use:

The SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER must be removed from its packaging and handled in aseptic conditions.

Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact.

The SINGLE-USE 0.3 PUNCTUM DILATOR AND PLUG INSERTER is a single-use product and must not be resterilized. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

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