REF. 2138 296 79 F

FR	NOTICE D'UTILISATION - SET DE GABARITS (ENUCLEATION-EVISCERATION) A USAGE UNIQUE	3
EN	INSTRUCTION FOR USE - SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION)	4
DE	GEBRAUCHSANWEISUNG - SCHABLONENSET (ENUKLEATION-EVISCERATION) ZUM EINMALGEBRAUCH	5
IT	ISTRUZIONI PER L'USO - SET DI CALIBRI (ENUCLEAZIONE-EVISCERAZIONE) USA E GETTA	6
ES	INSTRUCCIONES DE USO - SET DE PLANTILLAS (ENUCLEACIÓN-EVISCERACIÓN) PARA USO ÚNICO	7
PT	INSTRUÇÕES DE USO - SET DE GABARITOS (ENUCLEAÇÃO-EVISCERAÇÃO) PARA USO ÚNICO	8
NL	GEBRUIKSAANWIJZING - KIT VAN KALIBERS (ENUCLEATIE-EXENTERATIE) VOOR EENMALIG GEBRUIK	9
DA	BRUGSVEJLEDNING - STØRRELSESMÅLER TIL ENGANGSBRUG (ENUKLEATION/EVISCERATION)	10
FI	KÄYTTÖOHJEET - KERTAKÄYTTÖINEN TULKKISETTI (ENUNKLEAATIO-EVISKERAATIO)	11
SE	BRUKSANVISNING - SCHABLONSET (ENUKLEATION-EVISCERATION) FÖR ENGÅNGSBRUK	12
TR	KULLANIM - TEK KULLANIMLIK AYIRICI SETİ (ENÜKLEASYON/EVİSERASYON)	13
ZH	使用说明 - 一次性筛分器套件 (眼球摘除术 - 眼内容物剔除)	14
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΣΕΤ ΜΕΤΡΗΤΩΝ ΜΙΑΣ ΧΡΗΣΗΣ (ΕΞΟΡΥΞΗ-ΕΞΕΝΤΕΡΩΣΗ)	15
PL	INSTRUKCJA UŻYTKOWANIA - JEDNORAZOWY ZESTAW PRZYMIARÓW (ENUKLEACJA-EWISCERACJA)	16
SL	NAVODILA ZA UPORABO - KOMPLET MERILNIKOV ZA ENKRATNO UPORABO (ENUKLEACIJA-EVISCERACIJA)	17
SK	NÁVOD NA POUŽITIE - SET DE GABARITS (ENUKLEÁCIA-VIZCERÁCIA) A USAGE UNIQUE	18

Symbols 19

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
EN

Instructions for use
Date of revision of the instructions for use : 03/2016

SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION)

Description :

The SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION) includes 4 sizes, 14, 16, 18 and 20 mm in diameter, made of polyoxymethylene.

Packaging :

The SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION) is packaged in a double packaging and supplied sterile. It is sterilized in its final packaging using ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green.

Indications :

The SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION) is designed to determine the proper orbital implant diameter for enucleation or evisceration.

Recommendations for use :

- The diameter choice of an ORBITAL IMPLANT is an important parameter for the final result.
- Given the retractile effect of orbital tissues, and to avoid dehiscence around the scar, the ORBITAL IMPLANT diameter has to be adapted to the patient's anatomy.
- The SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION) is designed to determine the proper orbital implant diameter for enucleation or evisceration.
- The proper diameter should allow the wound to be closed without tension. The choice of the diameter must take into account the thickness of any material used to wrap the implant.

Precautions for use :

The SINGLE-USE SIZERS SET (ENUCLEATION-EVISCERATION) must be removed from its packaging and handled in aseptic conditions.

Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact.

The SIZERS SET (ENUCLEATION-EVISCERATION) is a SINGLE-USE product and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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