REF. 2138 315 22 A.1

FR	NOTICE D'UTILISATION - RITLENG[®]+	2
EN	INSTRUCTION FOR USE - RITLENG[®]+	4
DE	GEBRAUCHSANWEISUNG - RITLENG[®]+	6
IT	ISTRUZIONI PER L'USO - RITLENG[®]+	8
ES	INSTRUCCIONES DE USO - RITLENG[®]+	10
PT	INSTRUÇÕES DE USO - RITLENG[®]+	12
NL	GEBRUIKSAANWIJZING - RITLENG[®]+	14
DA	BRUGSVEJLEDNING - RITLENG[®]+	16
FI	KÄYTTÖOHJEET - RITLENG[®]+	18
SE	BRUKSANVISNING - RITLENG[®]+	20
TR	KULLANIM - RITLENG[®]+	22
ZH	使用说明 - RITLENG[®]+	24
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - RITLENG[®]+	26
PL	INSTRUKCJA UŻYTKOWANIA – RITLENG[®]+	28
SL	NAVODILA ZA UPORABO – RITLENG[®]+	30
SK	NÁVOD NA POUŽITIE - RITLENG[®]+	32

Figure **34**

Symbols **35**

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RITLENG®+**Description and packaging:**

The RITLENG®+ (S1.1455, S1.1456) is an autostable nasal bicanalicular intubation for the lacrimal ducts. It consists of a silicone tube with two larger diameter portions, which is connected at each extremity with a PEEK thread guide. Each thread guide has a wide diameter part and a narrow diameter part. The silicone tube may be coated with polyvinylpyrrolidone (PVP) to improve its wettability (RITLENG®+ PVP S1.1495, S1.1496). The central body of the intubation displays a central marking intended for the intubation placement between the lacrimal puncta.

The RITLENG®+ intubations are supplied sterile. They are sterilized by ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. The RITLENG®+ intubations are sterilized in their final double packaging to make them easier to handle in aseptic conditions.

Indications:

Nasal bicanalicular intubation is indicated in the treatment of epiphora in patients of 12 months and older in cases of:

- Canalicular pathologies (stenoses, obstructions, lacerations),
- Congenital duct obstruction (Hasner valve stenosis),
- Dacryocystorhinostomy (conventional or laser).

Mode of action:

The silicone tube acts as a conformer and enables drainage of tears by capillarity. In cases of canalicular lacerations, the silicone tube guides wound healing and reduces post-traumatic shrinkage. The advantages of the RITLENG®+ intubation over conventional bicanalicular probes are the following :

Removing a PEEK thread is less traumatic to the nasal mucosa, less haemorrhagic and less painful under local anesthesia, than removing a metal guide. Locating a PEEK buckle that lies wide and flat in the nose is much easier than locating a metal guide that comes out in the inferior meatus, often difficult to reach.

Warnings and precautions for use:

As with all other forms of intubation, a visual monitoring of the nasal cavity is recommended with the help of intranasal illumination or a nasal endoscope (0° or 30°) which makes it easier to locate the thread guide, particularly if it takes a posterior direction.

It is recommended not to remove the PEEK when the guide probe is in the nasolacrimal duct as it can get hooked in the groove of the probe becomes fragile and tears.

Contra-indications:

The RITLENG®+ intubations shall not be used in the following cases:

- Lacrimal sac tumor,
- Dacryolith,
- Other causes of tearing without a lacrimal origin, including but not limited to, tear hypersecretion (for example, dry eyes, blepharitis, allergic conjunctivitis, corneal ulceration, etc); lacrimal pump failure (for example facial paralysis).

Adverse side effects:

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of a RITLENG®+ intubation are similar to those associated with conventional nasal bicanalicular intubations, the most frequent being :

Intraoperative complications :

- false passages
- separation of the tube and the PEEK thread
- difficulty in locating the thread in the nose (hence the importance of endoscopy).

Postoperative complications (rare):

- conjunctival or nasal pruritus
- irritation in the internal angle of the eye or nasal irritation
- reversible narrowing of the palpebral fissure
- exteriorization or loss of the probe
- stricturotomy
- synechia of the nasal mucosa
- canaliculitis
- induced mucocele

A coloration of the probes is possible with the use of eye drops. It is not a risk regarding the tolerance of the probe.

Unexpected side effects and complications related to the RITLENG®+ intubations must be reported to FCI.

Instructions for use:Insertion

Using the RITLENG®+ intubation requires use of the RITLENG® PROBE (S1-1460) which is specially designed for this purpose. The intubation technique presents the following particularities:

Probing of the first canaliculus towards the nasolacrimal duct, up to the floor of the nasal fossa.

After this probing with the RITLENG® PROBE and guiding its plateau to direct the inferior opening towards the front of the nose, the PEEK thread guide is introduced in the RITLENG® PROBE so as to obtain a wide loop which will lie flat on the nasal fossa.

It is strongly recommended that a nasal endoscope and fine endonasal forceps be used for easy recovery of the PEEK thread.

In the absence of endoscopy, take a small blunt hook, seek contact with the inferior extremity of the metal probe and rotate the guide probe towards the rear. Then gently retract the probe and as soon as contact with the metal is lost, pull gently on the hook and remove the thread.

The PEEK thread (wide section) is removed from the nose and the metal probe is extracted from the lacrimal ducts.

The RITLENG® PROBE is then removed from the intubation by pending the narrow section of the PEEK thread guide out through the slit.

The same procedure is followed to probe the second canaliculus to achieve a nasal bicanalicular intubation, with an internal angle loop that is neither too loose nor too tight. The two exposed guide threads allow the setting of the intubation and its centering. The surplus of silicone tube is cut and removed.

Verify the correct positioning of the probe by asking the patient to blink several times.

It is recommended that the patient be shown the interpalpebral loop after the procedure.

Removal:

Grasp the visible loop with a blunt forceps.

With a slight but firm pressure, remove an end of the probe by withdrawing it from the lacrimal punctum.

Grasp the second end of the probe and remove it firmly.

The RITLENG®+ intubation must be removed from its packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. The RITLENG®+ intubations are single-use products and must not be resterilized; reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. Do not use the product if the indicator is not green as it might endanger the patient. They should be stored at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 4 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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