**REF. 2138 309 41 C**

<b>FR</b>	<b>NOTICE D'UTILISATION - SONDES DE CRAWFORD</b>	<b>2</b>
<b>EN</b>	<b>INSTRUCTION FOR USE - CRAWFORD INTUBATION</b>	<b>3</b>
<b>DE</b>	<b>GEBRAUCHSANWEISUNG - CRAWFORD-SONDE</b>	<b>4</b>
<b>IT</b>	<b>ISTRUZIONI PER L'USO - SONDE DI CRAWFORD</b>	<b>5</b>
<b>ES</b>	<b>INSTRUCCIONES DE USO - SONDAS DE CRAWFORD</b>	<b>6</b>
<b>PT</b>	<b>INSTRUÇÕES DE USO - SONDAS DE CRAWFORD</b>	<b>7</b>
<b>NL</b>	<b>GEBRUIKSAANWIJZING - SONDES VAN CRAWFORD</b>	<b>8</b>
<b>DA</b>	<b>BRUGSVEJLEDNING - CRAWFORD-SONDER</b>	<b>9</b>
<b>FI</b>	<b>KÄYTTÖOHJEET - CRAWFORD-SONDIT</b>	<b>10</b>
<b>SE</b>	<b>BRUKSANVISNING - CRAWFORD-SOND</b>	<b>11</b>
<b>TR</b>	<b>KULLANIM - CRAWFORD PROBLARI</b>	<b>12</b>
<b>ZH</b>	<b>使用说明 - CRAWFORD 探针</b>	<b>13</b>
<b>EL</b>	<b>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - ΔΙΑΣΩΛΗΝΩΣΗ CRAWFORD</b>	<b>15</b>
<b>PL</b>	<b>INSTRUKCJA UŻYTKOWANIA ZESTAW DO INTUBACJI CRAWFORDA</b>	<b>17</b>
<b>SL</b>	<b>NAVODILA ZA UPORABO – INTUBACIJA CRAWFORD</b>	<b>18</b>
<b>SK</b>	<b>NÁVOD NA POUŽITIE - CRAWFORDOVA INTUBÁCIA</b>	<b>19</b>

**Figures ..... 20**

**Symbols ..... 21**

**FCI S.A.S.**

20-22, rue Louis Armand 75015 PARIS  
- FRANCE

Tél : + 33 (0) 1 53 98 98 98

Fax : + 33 (0) 1 53 98 98 99

Email : [fci@fci.fr](mailto:fci@fci.fr) Web :  
[www.fci.fr](http://www.fci.fr)

Distributed in U.S.A. by :  
FCI Ophthalmics Inc. 20 Winter  
Street  
Pembroke, MA 02359  
Phone: 800-932-4202  
Fax: 781-826-9062

Email : [info@fci-ophthalmics.com](mailto:info@fci-ophthalmics.com) Web :  
[www.fci-ophthalmics.com](http://www.fci-ophthalmics.com)

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## CRAWFORD INTUBATION

### Description :

The CRAWFORD INTUBATION stents are composed of a silicone tube that is attached at each end to a flexible metallic probe with a diameter of 0.4 mm ending with an olive tip measuring 0.9 mm in diameter.

### Packaging :

The CRAWFORD INTUBATION stents are supplied sterile. They are sterilized by ethylene oxide. A green indicator shows that the product followed a sterilization cycle validated by FCI. Do not use the product if the indicator is not green. The CRAWFORD INTUBATION stents are sterilized in their final double packaging to make them easier to handle in aseptic conditions.

### Indications :

The CRAWFORD INTUBATION stents are indicated in treatments of epiphora in patients 12 months of age and older. Indications for nasolacrimal intubation performed with the CRAWFORD INTUBATION stents are:

- Canalicular pathologies (stenosis, obstruction, lacerations),
- During dacryocystorhinostomy procedures (conventional or laser),
- Imperforation of the nasolacrimal duct in the infants.

### Mode of action :

Nasal bicanalicular intubations are used for intubation of the lacrimal ducts. The silicone tube acts as a conformer and enables drainage of tears by capillarity. In cases of canalicular lacerations, the silicone tube guides wound healing and prevents the onset of synechia.

### Warning and contraindications :

All precautions required for a surgical operation are to be followed. It is recommended to present to the patient and/or to their families, the interpalpebral loop after insertion. It is recommended to advise the patient and/ or their family that the probe must not be accidentally cut close to the eyelids.

The CRAWFORD INTUBATION stents should not be used in the following cases:

- Lacrimal sac tumor,
- Dacryolith

In the case of complications or unusual symptoms it is recommended to quickly contact an ophthalmologist.

### Adverse side effects :

As in any type of surgery, there are risks linked to the material or to developments of the initial pathology. Potential complications associated with the implantation of the CRAWFORD INTUBATION stents include, but are not limited to the following :

#### Intraoperative complications :

- |   |                      |
|---|----------------------|
| - false passages                                  | - bleeding           |
| - separation of the silicone from the steel guide | - canalicular trauma |

#### Post-operative complications :

- |   |                                |
|---|--------------------------------|
| - conjunctival or nasal pruritus                | - synechia of the nasal mucosa |
| - nasal or caruncular irritation                | - canaliculitis                |
| - reversible shrinkage of the palpebral fissure | - induced mucocele             |
| - exteriorization or loss of the stent          | - bleeding                     |
| - stricturotomy                                 |                                |

Unexpected side effects and complications related to the CRAWFORD INTUBATION stents must be reported to FCI.

### Recommendations for use :

#### *Insertion*

Each punctum should be dilated. The lacrimal system is probed to open any blockages by using a standard lacrimal probe. The first probe of the Crawford Bicanaliculus Intubation stent is passed through the upper punctum and across the upper canaliculus and then oriented down through the lacrimal system into the nose approximately 4 cm. Proper positioning of the stent probe is confirmed by achieving metal-to-metal contact with the Crawford Hook (S1.1275) inserted through the nostril. The stent will be located lateral to the inferior turbinate in the inferior meatus of the nose. The hook engages the stent probe at its olive tip end and pulls the stent out of the nose. The second stent probe is passed down the inferior canaliculus and out of the nose in a similar fashion. The two exposed stent probes are cut and removed and the silicone tubes are tied to each other and securely knotted. The silicone tubes may be further secured by a fine suture of the surgeon's choice. The knotted silicone tubes are then tucked up into the nose.

#### *Removal*

Removal of the silicone tubes is accomplished by locating the silicone tube in the eye between the upper and lower puncta and pulling the tubing upward. The silicone tube is then cut and pulled out completely from the ducts. If the knot is too big and does not pass through the canaliculus, the silicone tubes must be removed through the nose. In older patients, it may be possible to expel the knot by having the patient blow his nose. In small children, location of the knot and silicone tube removal may need to be performed under general anesthesia.

The silicone - steel guide junction is a fragile area. It is recommended that the steel guide be held in front of the junction while handling the probe.

The use of this probe may require the use of a Crawford hook (S1.1275) that has been specifically designed for the lancet.

### Precautions for use :

CRAWFORD INTUBATION stents must be removed from their double packaging and handled in aseptic conditions. Before use, the individual pack preserving the product's sterility should be checked to make sure it is intact. CRAWFORD INTUBATION stents are single-use products and must not be resterilized. Reuse and/or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety. They should be stored at room temperature and must not be used after the expiration date shown on the package.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

An implant card is available in the product's packaging. It is intended to be provided to the patient and allows traceability of the implanted product.

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