REF. 2138 312 91 C

FR	NOTICE D'UTILISATION - RETILOCK®	2
EN	INSTRUCTION FOR USE - RETILOCK®	3
DE	GEBRAUCHSANWEISUNG - RETILOCK®	4
IT	ISTRUZIONI PER L'USO - RETILOCK®	5
ES	INSTRUCCIONES DE USO - RETILOCK®	6
PT	INSTRUÇÕES DE USO - RETILOCK®	7
NL	GEBRUIKSAANWIJZING - RETILOCK®	8
DA	BRUGSVEJLEDNING - RETILOCK®	9
FI	KÄYTTÖOHJEET - RETILOCK®	10
SE	BRUKSANVISNING - RETILOCK®	11
TR	KULLANIM - RETILOCK®	12
ZH	使用说明 - RETILOCK®	13
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ - RETILOCK®	14
PL	INSTRUKCJA UŻYTKOWANIA – RETILOCK®	15
SL	NAVODILA ZA UPORABO – RETILOCK®	16
SK	NÁVOD NA POUŽITIE - RETILOCK®	17

Symbols 18

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EN

Instructions for use
Date of revision of the instructions for use : 05/2016

RETILOCK®

Description :

The RETILOCK® TROCARS kits are subdivided into several reference numbers (box of 5):

S9.7000.23 : RETILOCK® 23G NON VALVED SELF-RETAINING TROCARS

S9.7100.23 : RETILOCK® 23G VALVED SELF-RETAINING TROCARS

S9.7000.25 : RETILOCK® 25G NON VALVED SELF-RETAINING TROCARS

S9.7100.25 : RETILOCK® 25G VALVED SELF-RETAINING TROCARS

They are composed of the following elements:

- Three cannulas equipped with an annular bulge, thus reducing the risk of accidental removal of each cannula. They are mounted on the blade for "One step" use. The knives are equipped with a scleral marker located at the rear of the handle in order to mark the distance at 3 or 4mm from the limbus.
- An infusion line can be connected to the cannulas to supply physiological saline solution or BSS, in order to maintain the tonus of the eye during surgery. The infusion line is compatible with the injection of silicone oil 1000 or 5000cst, which is necessary for the internal tamponade.
- The valved version comprises a valve mounted on the cannula, and the non-valved version comprises three plugs that can be adjusted on the cannula. These elements ensure that the cannula remains sealed when no instrument is inserted.

Presentation :

- The RETILOCK® TROCARS kits are supplied in sterile double blisters to facilitate handling under aseptic conditions.
- They are sterilized by ethylene oxide in their final packaging.
- A green indicator indicates that the product has undergone a sterilization cycle.
- Do not use the product if the indicator is not green.

Precautions for use:

- The devices must be removed from their blisters and handled under aseptic conditions.
- Check the integrity of the individual sterility protector before use.
- The RETILOCK® TROCARS kits are single-use devices that must not be resterilized.
- Reuse and/or resterilization may affect the performance of the device and cause a risk of compromising the safety and health of the patient.
- The devices must be stored at room temperature and must not be used after the expiry date stated on the packaging.

Indications :

The RETILOCK® TROCARS kits are intended for transconjunctival sclerotomy in order to gain access to the posterior pole.

Warnings:

- In the non-valved version, the cannulas are sealed by means of the plugs. Sealing is obtained by the weight of the plug itself, simply placing it in the cannula. It is not necessary to press on the plugs when they are inserted into the cannulas; this could block them. If the plug is blocked, slide a forceps between the plug head and the cannula head, and tighten the clamp to detach it from the cannula.
- During the injection of silicone oil using the infusion line included in the kit, the operating pressure must not exceed 5 bar (72.5 psi).
- It is the responsibility of the user to ensure compatibility between the internal lumen of cannulas of RETILOCK® and the size of the instruments.
- Only the infusion lines 23G and 25G (supplied with the kit RETILOCK® TROCARS) are respectively compatible with 23G and 25G cannulas of RETILOCK® TROCARS. The use of an infusion line from another manufacturer may compromise performance RETILOCK® TROCARS kits and patient safety.

Advice for use:

- It is recommended to use the scleral marker located at the rear of the knife handle in order to mark the limbus at 3 to 4 mm before making the incision (optional).
- At the end of the procedure, the surgeon may decide at his discretion whether the incisions should be sutured.

The guaranteed maximum residual rate of ethylene oxide at the time the device was placed on the market was 5.3 µg per device. The residual quantity of ethylene oxide corresponds to the maximum quantity of residue that may leach over the lifetime of use.

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